



CINCAP GmbH
Switzerland



IFIIP GmbH
Switzerland

“Right, First Time” Concept & Workflow based on F-CAD (Formulation Computer-Aided Design)

Prof. Dr. Hans Leuenberger

**IFIIP GmbH, Institute for innovation in industrial pharmacy
Kreuzackerweg 12, CH-4148 Pfeffingen, Switzerland**

www.ifiip.ch

Dr. Maxim Puchkov

**CINCAP GmbH, Center for innovation in computer-aided pharmaceuticals
Kreuzackerweg 12, CH - 4148 Pfeffingen, Switzerland**

www.cincap.ch

Dr. Bernhard Schneider

RPD Tool AG, Rothausstrasse 61, CH-4132 Muttenz, Switzerland

www.rpdtool.com



CINCAP GmbH
Switzerland




IFIP GmbH
Switzerland

Formulation Computer-Aided Design (F-CAD): Key for “Right, First Time” Workflow to speed -up Time to Market

The drug development process

Graph shows declining product approval despite increasing R&D investment. R&D investment figures only include company-financed R&D. Source: *PhRMA member survey 2011 and Nature Drug Discovery*¹



¹ Mullard, A, *Nature Reviews Drug Discovery* 10, 82-85 (February 2011)




CINCAP GmbH
Switzerland



IFIIP GmbH
Switzerland

Attrition rate of drugs tested





CINCAP GmbH
Switzerland



IFIP GmbH
Switzerland

*Formulation Computer-Aided Design (F-CAD):
Key for “Right, First Time” Workflow
to speed -up Time to Market*

Sharp rise in the attrition rate between 1990 and 2004*):

Drug companies are removing in Clinical Phase III
more compounds from the pipeline at all levels of testing!

Different reasons: 1) Lack of Efficacy 2) Safety 3) Commercial/
Financial 4) Not disclosed

Lack of Quality by Design at all levels?

*) S. Schreder, SWISS PHARMA 34 (2012) p.19



CINCAP GmbH
Switzerland



IFIP GmbH
Switzerland

*Formulation Computer-Aided Design (F-CAD):
Key for “Right, First Time” Workflow
to speed -up Time to Market*

Sharp rise in the attrition rate between 1990 and 2004*):

Drug companies are removing in Clinical Phase III
more compounds from the pipeline at all levels of testing!

Different reasons: 1) Lack of Efficacy 2) Safety 3) Commercial/
Financial 4) Not disclosed

Lack of Quality by Design at all levels?

*) S. Schreder, SWISS PHARMA 34 (2012) p.19




CINCAP GmbH
Switzerland



IFIP GmbH
Switzerland

FDA Whitepaper March 2004

Three Dimensions of the Critical Path






CINCAP GmbH
Switzerland



IFIP GmbH
Switzerland

Quality by Design

Product and Process Quality Knowledge: Science-Risk Based cGMP's





CINCAP GmbH
Switzerland



IFIP GmbH
Switzerland

Conventional workflow for clinical trials

Early Phase: Service dosage form (in general capsules)



Change to tablets & Bioequivalence Testing



Scale-Up Exercise



Mass-production of final marketed form (Two-sigma quality)



CINCAP GmbH
Switzerland



IFIP GmbH
Switzerland

Right-First-Time Workflow

Market ready tablet dosage form (instead of service form)



Small-scale Production



Scale-Up Exercise (computer-assisted)



Mass-production of final marketed form (six-sigma quality)



CINCAP GmbH
Switzerland



IFIP GmbH
Switzerland

*Formulation Computer-Aided Design (F-CAD):
Key for “Right, First Time” Workflow
to speed -up Time to Market*

Regulatory Issues:

ICH Q 8 (R2)

***Design space exploration with F-CAD
eliminates risk for
undesired time extension***




CINCAP GmbH
Switzerland




IFIP GmbH
Switzerland

Conventional versus Right First Time (RFT) Workflow

Conventional workflow



RFT workflow





CINCAP GmbH
Switzerland



IFIP GmbH
Switzerland

F-CAD

Computer-Aided Formulation Design



F-CAD is an ultimate set of modeling and computational tools to assist in formulation design of pharmaceutical solid dosage forms with the goal to save money by replacing lab work with “in-silico” experiments.



Anticipated market ready tablet formulation

Steps:

1)

- > Selection of the drug substance:
- > for a comprehensive test-> check all classes of BCS
- > Physico-chemical properties of the drug substance needed

- > Choice of chemical/physical compatible excipients
- > Choice of the strength of the tablet formulation, e.g. 5 strengths
- > Choice of the dissolution profile ->immediate release and/or -> controlled release

2)

- > Application of F-CAD
- > In-silico evaluation of the Design Space according to ICH Q8 (R2): Sensitivity Report (excipients, process variables) on dissolution profile

3)

- >Choice of the best in-silico dosage form (different strengths) for Phase 1 Clinical Study for laboratory validation and clinical manufacturing.
- > Normal chem. stability program (incl. Dissolution rate, disintegr. Time)
- > Accelerated Stability Test Program for Formulations: 1. results after 1-2 weeks.







CINCAP GmbH
Switzerland



IFIP GmbH
Switzerland

F-CAD Tablet Designer



3D Printer, Tooling Manufacturer,
etc. OR export to F-CAD PAC




CINCAP GmbH
Switzerland



IFIP GmbH
Switzerland

Custom Shapes with F-CAD Tablet Designer







CINCAP GmbH
Switzerland



IFIP GmbH
Switzerland

F-CAD PAC – Particle Arrangement and Compaction






CINCAP GmbH
Switzerland



IFIP GmbH
Switzerland

Resulting compact



**First Principle,
i.e. Cellular Automata
approach for calculating
dissolution rate profile**




CINCAP GmbH
Switzerland



IFIIP GmbH
Switzerland

Experimental vs. in silico dissolution profiles of different formulations with caffeine: Experimental * and in-silico formulation standard error due to random arrangement of drug particles in the tablet






CINCAP GmbH
Switzerland




IFIP GmbH
Switzerland

F-CAD – *in silico* Dissolution Profiles: 1) Reproducibility of details (see zoom) 2) Effect of tablet shape




© Heath Scientific



↑ 1) Initial dissolution detail

← 2) Effect of tablet shape





Percolation Theory and F-CAD

F-CAD is capable to detect Percolation thresholds in a tablet formulation, which can be the source of the variability of a tablet property such as the disintegration time, see PhD Thesis Go Kimura, eLink: http://edoc.unibas.ch/diss/DissB_9886




Fig. 68 of PhD thesis: Experimental disintegration time (●) and F-CAD values (△) as a function of the renormalized MS (= Maize starch as disintegrant) concentration.

!



CINCAP GmbH
Switzerland



IFIIP GmbH
Switzerland

F-CAD and Quality by Design

- » **F-CAD enables and facilitates to comply with the requirement of ICH Q8 (R2) at affordable costs!**
- » Screening for robust formulation
- » Setting up acceptance criteria for raw materials
- » Analyse scale-up/scale-down issues and prevent problems



CINCAP GmbH
Switzerland



IFIIP GmbH
Switzerland

Summary: Goals of F-CAD

- » Superior quality of formulations than with existing approach
- » Possibility to quantify the robustness of the formulation
- » Possibility to define specifications based on science
- » Reduction of time to market
- » Boosting formulation and process technology understanding
- » F-CAD is able to detect percolation thresholds within the design space exploration according to ICH Q8 (R2)
- » Computer aided design of formulations similar to aircraft design
- » Savings comparable to the savings of the aircraft industry



CINCAP GmbH
Switzerland



IFIP GmbH
Switzerland

Thank you for your attention!

Audience Q&A